

iStat Operator Competency

| |
|-------------------------|
| Name: |
| Employee Number: |
| Position: |
| Date: |
| Trainer: |

Operators to sign groups where indicated once completed

| | |
|---|-------------------|
| 1. Cartridges and Sample Requirements | User signature |
| <ul style="list-style-type: none"> Understands storage requirements for cartridges Checks active expiry of cartridge Aware of sample type requirements for each cartridge Ability to perform electronic simulator test | |
| 2. Sample Testing | User signature |
| <ul style="list-style-type: none"> Turn on the analyser Scan User Login Scan encounter number from patient Bradma label Scan the correct cartridge lot number Handles test cartridge correctly Mixes sample adequately Applies sample correctly to cartridge Inserts cartridge and waits for results Disposal of cartridge Prints result Records results appropriately Can retrieve a patient result Basic troubleshooting knowledge Change/replace battery Change printer paper | |
| 3. Available resources | User signature |
| <ul style="list-style-type: none"> Reviewed user manual and quick guides | |
| 4. Quality Control - Lead Users Only | User signature |
| <ul style="list-style-type: none"> Know when to run QC – Monthly and on new Lots/Deliveries Running, recording, expiry & storage QC Understanding and running of EQA | |
| Signed by Trainer (POCT QM/Educator/Lead User) | Trainer signature |