

Clinitek Urinalysis Operator Competency Assessment

Name:
Employee Number:
Position:
Date:
Trainer:

Operators to sign groups where indicated once completed

1. Sample Collection & Consumables	User signature
<ul style="list-style-type: none"> Aware of Health and Safety requirements Aware of acceptable and unacceptable samples Aware that only Siemens Multistix can be used on instrument 	
2. Sample Testing	User signature
<ul style="list-style-type: none"> Turn on/off the analyser Enter User, Patient, Strip details Perform patient test correctly and timely Understand the importance of timing for accurate results Discard strip and wipe test tray Record result in patient notes Can retrieve a patient result How to change paper/stickers Knowledge of QC requirements 	
3. Available Resources	User signature
<ul style="list-style-type: none"> Reviewed user manual 	
4. Quality Control - Lead Users Only	User signature
<ul style="list-style-type: none"> Know when to run QC – Weekly and on new Strip Perform QC correctly Understanding of EQA Know who & when to contact for assistance 	
Signed by Trainer (POCT QM/Educator/Lead User)	Trainer signature